

PRODUCT SPECIFICATION

REF 6036

SALVEQUICK PLASTIC PLASTERS (ADHESIVE BANDAGES)



1 GENERAL PRODUCT DESCRIPTION

Description

Adhesive bandages (plasters) of two sizes, medium and senior, packed in a cardboard refill package, intended for the Cederroth (Salvequick) dispenser for adhesive bandages.

The plasters can easily be pulled out with just one hand and applied on the wound.

The adhesive bandage consists of a backing, evenly spread with an adhesive coating. An absorbent pad is fixed centrally to the dressing. The adhesive is covered with a release-coated foil that consists of two parts partly overlapping each other. When the adhesive bandage is taken out from the package, the release-coated plastic foil remains behind. One half of the dressing can now be applied on the skin using one hand only. The rest of the coated paper can then be removed and the adhesive bandage is ready for use.



Storage information

Dry storage in room temperature. For maximum protection, keep in original package until use.

2 DETAILED PRODUCT INFORMATION

Content per consumer unit

<u>Product/Material</u>	<u>Quantity (pcs)</u>	<u>Dimensions (mm)</u>	
Refill packages	6	30 x 88 x 100	
<u>Per refill package:</u>		<u>Dressing</u>	<u>Wound Pad</u>
Medium plasters	27	72 x 19	24 x 12
Senior plasters	18	72 x 25	24 x 17

Material

Film/backing:	Beige perforated polyethylene foil
Adhesive:	Acrylate
Wound pad:	Viscose, polypropylene, polyethylene with polyethylene wound contact layer
Facing:	Siliconized paper and plastic
Protective paper:	Cold seal paper coated with latex-based adhesive
Wrapping:	Carton

Shelf life: 5 years

3 LOGISTIC INFORMATION

	CU	SKU	Pallet
GTIN:	7310610060367	7310806036213	7333460013475
Content:	6 refills	10 CU	10 SKU/layer x 5 layers
Size (w x l x h mm):	127 x 102 x 90,5	139 x 537 x 206	Pallet 800 x 1200 x 1174
Gross weight:	198 g	2216 g	136,1 kg
Packing material:	Carton 32 g (>80% recycled fiber)	Corrugated board 236 g (100% recycled fiber)	-

