

PRODUCT SPECIFICATION

REF 51030126 SALVEQUICK DETECTABLE FINGERTIPS/REGULAR



1 GENERAL PRODUCT DESCRIPTION

Description

Metal Detectable:

Detectable, waterproof, adhesive bandages (plasters), of two sizes, medium and fingertip, packed in a cardboard refill intended for the Cederroth (Salvequick) dispenser for adhesive bandages.

The Salvequick Detectable bandages are detectable by most common metal detectors and the blue color also makes them easy to detect visually, which makes these plasters suitable for restaurants and the food industry. The plasters can easily be pulled out with just one hand and applied on the wound.

The adhesive bandages consist of a backing, evenly spread with an adhesive coating and with an absorbent pad fixed centrally to the dressing. When the adhesive bandage is pulled out from the package, the lower cover will remain behind. Thus, one half of the bandage can be applied with one hand and after removing the rest of the cover the plaster is ready for use.

Storage information

Dry storage in room temperature. For maximum protection, keep in original package until use.



2 DETAILED PRODUCT INFORMATION

Content per consumer unit

<u>Product/Material</u>	<u>Quantity (pcs)</u>	<u>Dimensions (mm)</u>	
Refill packages	6	30 x 88 x 100	
<u>Per refill package:</u>		<u>Dressing</u>	<u>Wound Pad</u>
Medium plasters	15	72 x 19	24 x 12
Fingertip	15	66 x 39	28 x 12

Material

Film/backing:	Blue, soft and flexible polyurethane/polystyrene foil
Adhesive:	Acrylate
Wound pad:	Laminate of viscose, polypropylene, polyethylene with polyethylene wound contact layer and a metal detectable aluminium foil
Facing:	Siliconized paper and plastic
Protective paper:	Cold seal paper coated with latex-based adhesive
Wrapping:	Carton

Shelf life: 5 years

3 LOGISTIC INFORMATION

	CU	SKU	Pallet
GTIN:	7310610019310	7310800019397	7333460013581
Content:	6 refills	10 CU	10 SKU/layer x 5 layers
Size (w x l x h mm):	127 x 102 x 90,5	139 x 537 x 206	Pallet 800 x 1200 x 1174
Gross weight:	173 g	1966 g	123,6 kg
Packing material:	Carton 32 g (>80% recycled fiber)	Corrugated board 236 g (100 % recycled fiber)	-

