



Cederroth First Aid

Blue detectable plasters for the food industry

Detectable plasters are required when handling food

- ✓ Visually recognisable by the blue colour
- ✓ Metal and/or X-Ray detectable



SALVEQUICK X-Ray Detectable Plaster

New!

REF 510003320

Salvequick Blue X-Ray Detectable Plasters, specially designed for use in the food industry. The plasters are both X-ray and metal detectable.

- X-ray detectable
- Metal detectable
- Visually detectable
- Waterproof
- Easy and quick to apply
- Flexible and convenient
- Each refill contains 30 plasters (18 pcs 72x19 mm and 12 pcs 72x25 mm). 6 refills/box



SALVEQUICK X-Ray Detectable Plaster Dispenser

New!

REF 510003295

Salvequick X-Ray Detectable Plaster Dispenser with blue x-ray and metal detectable plasters is specially designed for use in the food industry. For quick and easy access to plasters. The plasters are visible and clean until needed.

- Easy to use
- Easy to refill
- Dust protected
- Contains X-Ray Detectable Plasters



Salvequick X-Ray Detectable Plaster Dispenser
REF 510003295
Dimensions: W 23 x H 12 x D 5,5 cm

2 (x 30) X-Ray Detectable Plaster (REF 510003320)
1 Refill key





SALVEQUICK

Blue Detectable Plaster & Plaster Dispenser

Salvequick Blue Detectable Plaster has been specifically designed for use in the food industry. The plasters are detectable by most metal detectors and the blue colour makes the plasters easier to see.

- Metal detectable
- Visually detectable
- Waterproof
- Easy and quick to apply
- Flexible and convenient
- Blue detectable plasters are available in two options



Salvequick Blue Detectable Plaster
REF 51030127

Each refill contains 35 plasters (21 plasters measuring 72 x 19 mm and 14 plasters measuring 72 x 25 mm). 6 refills/box.



Salvequick Blue Detectable Plaster Fingertip/Regular
REF 51030126

Each refill contains 30 plasters (15 plasters measuring 66 x 39 (22) mm and 15 plasters measuring 72 x 19 mm). 6 refills/box.



Salvequick Plaster Dispenser Blue Detectable
REF 51030130
Dimensions: W 23 x H 12 x D 5,5 cm

2 (x 35) Salvequick Blue Detectable Plasters (REF 51030127)
1 Refill key



In our Salvequick Dispenser you pull plasters downwards, which avoids you getting any dirt or blood on other plasters. The Salvequick Dispenser and Salvequick Plasters is a complete, easy-to-use system. The plasters are always available, ready to use. The plaster refills are locked in the dispenser to make sure they do not disappear. Empty refill packs are easy to remove using a special key.



Salvequick detectable products

New!



Salvequick X-Ray Detectable Plaster Dispenser

Content: 2 (x 30) X-Ray Detectable plaster. The plasters are both X-ray and metal detectable.
REF 510003295.

Suitable for all Cederroth
First Aid Dispensers and Stations



New!

Salvequick X-Ray Detectable Plaster

6 (x 30) X-Ray Detectable Plaster
REF 510003320.

Suitable for all Cederroth
First Aid Dispensers and Stations



Salvequick Blue Detectable Plasters

6 (x 35) Blue Detectable Plaster,
REF 51030127.
6 (x 30) Blue Detectable Plaster Fingertip/
Regular, REF 51030126.



Salvequick Blue Detectable Plaster Dispenser

Content: 2 (x 35) Blue Detectable Plaster. The plasters are metal detectable.
REF 51030130.



Wound Care Dispenser Blue

Content: 2 (x 35) Blue Detectable Plaster and 1 (x 30) Blue Detectable Plaster Fingertip/Regular. The plasters are metal detectable.
REF 51011009.



First Aid & Burn Station

Content: 1 (x 35) Blue Detectable Plaster and 1 (x 30) Blue Detectable Plaster Fingertip/Regular. The plasters are metal detectable.
REF 51011003.

Product		REF	CU / SKU	EAN	Launch
Salvequick X-Ray Detectable Plaster Dispenser inc. 2 (x 30) refills/box	 	510003295	5	7340213301765	1. June, 2024
Salvequick X-Ray Detectable Plaster 6 (x 30) refills/box	 	510003320	10	7340213301758	1. June, 2024
Salvequick Blue Detectable Plaster Dispenser inc. 2 (x 35) refills/box		51030130	20	7310610019693	
Salvequick Blue Detectable Plaster 6 (x 35) refills/box		51030127	10	7310610019327	
Salvequick Blue Detectable Fingertip/Regular 6 (x 30) refills/box		51030126	10	7310610019310	

Distributor:



Orkla Wound Care AB

(Svetsarvägen 15), Box 1336, SE-171 26 Solna, Sweden, Tel +46 (0)8 590 963 25. firstaid@cederroth.com www.cederroth.com