

PRODUCT SPECIFICATION

REF 510004108 SOFT FOAM BANDAGE 2 IN 1 DISPENSER



1 General product description

Description

The Cederroth Soft Foam Bandage 2-in-1 Dispenser is designed to accommodate two rolls in different sizes of blue Soft Foam Bandages, which are ideal for treating minor wounds.

The Soft Foam Bandage Blue is a self-adhering, adhesive-free plaster that remains in place even when exposed to water.

The dispenser can be securely mounted on the wall and features an integrated cutter for easy application.

The plastic used in the dispenser is made from 93% recycled materials, resulting in 90% less CO2 emissions during production compared to virgin plastic.



The dispenser contains CE-labelled Medical devices.

Storage information

Dry storage at room temperature. For maximum protection, keep it in its original package until use.

2 DETAILED PRODUCT INFORMATION

<u>Product/Material</u>	<u>Quantity</u>	<u>Dimensions</u>	<u>REF (refill)</u>
Dispenser	1	142x 138x 123 mm	-
Red key	1	-	-
Cederroth Soft Foam Bandage Blue 4,5m x 6cm Blue*	1	4,5m x 6cm	51011010
Cederroth Soft Foam Bandage Blue 4,5m x 3cm Blue*	1	4,5m x 3cm	51011022
Assembly instruction	1	-	-

**Soft Foam refills are also available in beige (REF 51011020, REF 51011018) and black (REF 51011021)*

Material

Dispenser:	Green parts:	>90% Recycled ABS
	Transparent lid and windows:	100% recycled PC
	Metal Parts:	
	Cutting board and rolling knife	Stainless steel
	Other plastic parts	
	Button	Yellow ABS
	Support disc	Natural color PA66
	Red key	Red Polyamide (PA)

Soft foam products: Refer to Product specifications for each product.

Shelf life: 5 years from production date



PRODUCT SPECIFICATION

REF 510004108 SOFT FOAM BANDAGE 2 IN 1 DISPENSER



3 LOGISTIC INFORMATION

	CU	SKU	Pallet
GTIN:	7340213302267	7333460000772	17333460001769
Content:	1 dispenser containing 2 pcs of blue soft foam bandages	6 CU	10 SKU/layer / 3 layers
Size (w x l x h mm):	145 x 140 x 138	165 x 445 x 300	Pallet 800 x 1200x 1045
Gross weight:	445 g	2 922 g	112,67 kg
Packing material:	Carton 73,71 g + 3 g Shrink-wrap	Carton 255g	-
	58% recycled fibers	82% recycled fibers	-

