

PRODUCT SPECIFICATION



REF 726000 CEDERROTH EYE & WOUND CLEANSING SPRAY

1 GENERAL PRODUCT DESCRIPTION

Description

A portable aerosol spray bottle filled with 150 mL physiological (0.9%) unbuffered sodium chloride solution, preservative free and sterile. The device provides a gentle spray douche that is suitable for use on both eyes and wounds. The actuator is protected by a plastic cap and a tamper seal. The nitrogen propellant is not in contact with the solution. The solution remains sterile until last drop. The Cederroth Eye & Wound Cleansing spray is intended to be used for rinsing and cleansing out dust and dirt from eyes and wounds.



Storage information

Store in temperatures between 15°C and 30°C.

2 DETAILED PRODUCT INFORMATION

Content per consumer unit

Material

Liquid	Sterile 0.9% Sodium Chloride (NaCl). Preservative free.
Can	Aluminum, 100% post-consumer recycled
Bag-on-valve	4 layer PP-OPA-ALU-PET
Actuator	White PP
Insert	Blue POM
Cover CAP	Transparent PP
Propellant	Nitrogen

Shelf life: 5 years from production date

3 LOGISTIC INFORMATION

	CU	SKU	Pallet
GTIN:	7310617260005	7310807260006	7312027260002
Content:	1 aerosol (150 ml)	12 CU	24 SKU/layer / 5 layers
Size (w x l x h mm):	53 x 53 x 162	168 x 221 x 191	Pallet 800 x 1200x 1105
Gross weight:	191 g	2,4 kg	317 kg
Packing material:	-	Carton 140 g	-
Transport information:	UN1950, non-flammable aerosols, class 2.2, Limited quantity LQ		

