

# PRODUCT SPECIFICATION

REF 390101 FIRST AID KIT MEDIUM



## 1 GENERAL PRODUCT DESCRIPTION

### Description

A First Aid kit containing first aid products to be used in emergency situations (bleeding, breathing, burns, sprains, circulatory shock). The strong and durable case has an adjustable handle making it easy to attach or carry around. All products are easily available and protected against dust, dirt and moisture inside the well-organised bag. There is also space for additional supplies.

The bright green and yellow coloured kit is clearly visible at a distance.

The Medical device products inside are CE-labelled.

### Storage information

Dry storage in room temperature. For maximum protection, keep the products in original package until use.



## 2 DETAILED PRODUCT INFORMATION

### Content per consumer unit

<u>Product/Material</u>	<u>Quantity</u>	<u>REF (refill)</u>
Cover	1	-
First Aid Instruction	1	51616402
Cederroth Breathing Mask	1	2596
Gloves	2 pairs	2596
Cederroth Safety Hand Cleanser	2	2596
Cederroth Emergency Blanket	1	1892
Scissors	1	270002
Cederroth 4 in 1 Bloodstopper, sterile	1	1910
Cederroth 4 in 1 Bloodstopper mini, sterile	1	1911
Cederroth Burn Gel dressing, sterile	1	901900
Salvequick Plasters	10	510005534
Salvequick Wound Cleansers, sterile	6	510005534
Salvequick Maxi Cover plaster (76 x 54 mm)	4	658024
Cederroth Elastic bandage	1	1882

### Material

Cover: PU coated EVA with 75D fabric inside. Inner pockets PE. Handle: Velcro (PP)+ Nylon  
Zipper: 300 D Nylon and plastic, puller made of zinc alloy. PES threads.

First Aid products: Refer to Product specifications for each product

Shelf life: 5 years from production date for each component

## 3 LOGISTIC INFORMATION

	CU	SKU	Pallet
<b>GTIN:</b>	7310613901018	7310803901019	7312023901015
<b>Content:</b>	1 kit	6 CU	9 SKU/layer x 4 layers
<b>Size (w x l x h mm):</b>	231 x 78 x 190	394 x 251x 240	Pallet 800 x 1200 x 1110
<b>Gross weight:</b>	576 g	3768 g	161 kg
<b>Packing material:</b>	Carton 20 g	Corrugated board 312 g	-



Version: 06  
Issued: 2026-05-11