

# PRODUCT SPECIFICATION

REF 390103 FIRST AID KIT X-LARGE



## 1 GENERAL PRODUCT DESCRIPTION

### Description

An extra large First Aid kit containing first aid products to be used in emergency situations (bleeding, breathing, burns, sprains, circulatory shock, minor eye injuries). The strong and durable case also has a Salvequick plaster dispenser on the outside allowing easy access to plasters without opening the case. Inside the well-organised case all products are easily available and protected against dust, dirt and moisture. The bright green and yellow coloured kit is clearly visible at a distance.

The Medical device products inside are CE-labelled.

**Accessory:** Wall bracket, REF 51000008



### Storage information

Dry storage in room temperature. For maximum protection, keep the products in original package until use.

## 2 DETAILED PRODUCT INFORMATION

### Content per consumer unit

<u>Product/Material</u>	<u>Qty</u>	<u>REF</u> <u>(refill)</u>	<u>Product/Material</u>	<u>Qty</u>	<u>REF</u> <u>(refill)</u>
First Aid Instruction	1	-	Metallized Dressing, 20x40cm, sterile	3	1893
Cederroth Breathing Mask	1	2596	Fixing Net	2	1893
Gloves, pairs	4	2596	Cederroth Burn Gel dressing, sterile	2	901900
Cederroth Safety Hand Cleanser	4	2596	Cederroth Soft Foam Bandage 6 x 40 cm	2	666150
Cederroth Emergency Blanket	1	1892	Cederroth 4 in 1 Bloodstopper, sterile	3	1910
Scissors	1	270002	Cederroth 4 in 1 Bloodstopper mini, sterile	4	1911
Tweezers	1	1875	Cederroth Elastic bandage 8cm x 4m	2	1882
Salvequick Plastic Plasters x 45	1	6036	Salvequick Cold Pack	1	219600
Salvequick Textile Plasters x 40		6444	Eye&Wound Cleansing Spray	1	726000
Salvequick Wound Cleansers, sterile	10	510005534			

### Material

Case: Case: HDPE, handle and plaster dispenser: PP, latches: POM, button lock and inside dividers: PS, Inside cover: acrylic glass.

First Aid products: Refer to Product specifications for each product

**Shelf life:** 5 years from production date for each component

## 3 LOGISTIC INFORMATION

	<b>CU</b>	<b>SKU</b>	<b>Pallet</b>
<b>GTIN:</b>	7310613901032	7310803901033	7312023901039
<b>Content:</b>	1 kit	2 CU	6 SKU/layer x 4 layers
<b>Size (w x l x h mm):</b>	420 x 118 x 300	430 x 305 x 249	Pallet 800 x 1200 x 1146
<b>Gross weight:</b>	3409 g	7150 g	197 kg
<b>Packing material:</b>	Carton 47 g	Corrugated board 332 g	-



Version: 04  
Issued: 2026-05-11