

PRODUCT SPECIFICATION

REF 546264 SALVEQUICK MED TEXTILE 5M

1 GENERAL PRODUCT DESCRIPTION

Description

Length plaster with an absorbing wound pad, intended to be used on minor wounds for protection and absorption of exudates.

The length plaster is elastic, breathable and has good adhesion.

CE-labelled Medical device.

Storage information

Dry storage in room temperature. For maximum protection, keep in original package until use.



2 DETAILED PRODUCT INFORMATION

Content per consumer unit

<u>Product/Material</u>	<u>Quantity</u>	<u>Dimensions (l x w mm)</u>	
		<u>Dressing</u>	<u>Wound pad</u>
Length plaster	5 m	500 x 60	500 x 21

Material

Film/backing:	Elastic, beige viscose/polyamide fabric with hydrophobic treatment
Adhesive:	Acrylate
Wound pad:	Viscose, polypropylene, polyethylene with polyethylene wound contact layer
Protective paper:	Siliconized paper

Shelf life: 5 years from production date

3 LOGISTIC INFORMATION

	CU	SKU	Pallet
GTIN:	7310615462647	7310805462648	7312025462644
Content:	5 m	4 CU	30 SKU/layer / 10 layers
Size (w x l x h mm):	98 x 75 x 98	196 x 150 x 98	Pallet 800 x 1200 x 1130
Gross weight:	107 g	432 g	155 kg
Packing material:	Carton 16 g	Plastic 4 g	-