

# PRODUCT SPECIFICATION

REF 1925 BREATHING MASK IN KEYRING POUCH



## 1 GENERAL PRODUCT DESCRIPTION

### Description

A breathing mask that consists of a one-way valve (green upper part, red lower part and a membrane) and a plastic sheet with graphic instructions printed on it. The Breathing mask is packed in a small pouch, which can be used as a key ring. The pouch is delivered in a plastic zip bag.



The device is intended to be used as an aid of assistance when giving mouth-to-mouth resuscitation to an injured person. The device is also intended to work as a barrier and limit the contact between the resuscitator and the injured person's blood, secrete, vomit and other body fluids.

CE-labelled.

### Storage information

Dry storage in room temperature. For maximum protection, keep in original package until use.

## 2 DETAILED PRODUCT INFORMATION

### Content per consumer unit

<u>Product/Material</u>	<u>Quantity</u>	<u>Dimensions (mm)</u>
Breathing mask (plastic sheet with valve)	1	27 ± 0,5 x 25 ± 0,5 Thickness 0,1 mm

### Material

Plastic sheet	PVC with phthalate (DEHP, DMP, BBP, DINP, DIDP, DNOP) contents < 0.003%
Valve, green upper part	Green ABS- plastic
Valve, red lower part	Red ABS-plastic
Valve, sealing	Red ABS-plastic
Valve, membrane	Red TPE (Thermoplastic polyester)
Plastic zip bag	Low density polyethylene

**Shelf life:** 5 years from production date

## 3 LOGISTIC INFORMATION

	CU	SKU	Pallet
<b>GTIN:</b>	7310611925016	7310801925000	7312020012196
<b>Content:</b>	1 piece	40 CU	14 SKU/layer x 6 layers
<b>Size (w x l x h mm):</b>	80 x 44 x 137	205 x 260 x 175	800 x 1200 x 1200
<b>Gross weight:</b>	20 g	976 g	107 kg
<b>Packing material:</b>	Plastic 2 g	Carton 155 g	-

