

## PRODUCT SPECIFICATION

**Product:** Refill for Cederroth First Aid Panel  
Panel

**Reference code:** 264000



### 1. General description

Refill products for Cederroth First Aid Panel (Ref 190900).

### 2. Intended use

For use in first aid situations.

### 3. Composition

#### *Consumer package*

<u>Ref#</u>	<u>Product/Material</u>	<u>Pcs</u>
1910	Cederroth 4-in-1	2
1911	Cederroth 4-in-1 mini	3
2596	Protection kit	1
6036CAP	Salvequick Plastic Adhesive bandages, refill	1
6444CAP	Salvequick Textile Adhesive bandages, refill	1
3227	Savett Wound Cleanser	1
51616410	First Aid Instruction	1
51616620	Key to the dispensers	1

#### *Transport package*

Number of consumer packages in transport package:	1
Dimensions (l x w x h) mm:	346 x 206 x 86
Weight (kg):	0,6

### 4. Labeling

#### *On transport package*

Cederroth address, telephone- and fax number, web-address

Product name, Reference code, Lot-number, EAN-code, Expiry date, Quantity

Symbols: LOT, USEBY

### 5. Storage conditions

The product should be stored in the original package for maximum protection until time of use.

Shelf life: 5 years from production date (Based on expiry date for components)

**Authorized by:** Teija Ålander  
*Quality Assurance Manager*

Date: 2008-11-25