



## PRODUCT SPECIFICATION

REF 6470 SALVEQUICK TEXTILE (FABRIC) XL  
PLASTERS (ADHESIVE BANDAGES)

### 1 GENERAL PRODUCT DESCRIPTION

#### Description

Adhesive bandages (plasters) of two sizes, small and big, packed in a cardboard refill package, intended for the Cederroth (Salvequick) dispenser for adhesive bandages.

The plasters can easily be pulled out with just one hand and applied on the wound.

The adhesive bandage consists of a backing, evenly spread with an adhesive coating. An absorbent pad is fixed centrally to the dressing. The adhesive is covered with a release-coated foil that consists of two parts partly overlapping each other. When the adhesive bandage is taken out from the package, the release-coated plastic foil remains behind. One half of the dressing can now be applied on the skin using one hand only. The rest of the coated paper can then be removed and the adhesive bandage is ready for use.

#### Storage information

Dry storage in room temperature. For maximum protection, keep in original package until use.



### 2 DETAILED PRODUCT INFORMATION

#### Content per consumer unit

<u>Product/Material</u>	<u>Quantity (pcs)</u>	<u>Dimensions (mm)</u>	
Refill packages	6	30 x 88 x 100	
<u>Per refill package:</u>			
		<u>Dressing</u>	<u>Wound Pad</u>
Small plasters	14	80 x 30	24 x 30
Big plasters	7	80 x 60	24 x 60

#### Material

Film/backing:	Elastic, beige viscose/polyamide fabric with hydrophobic treatment
Adhesive:	Acrylate
Wound pad:	Viscose, polypropylene, polyethylene with polyethylene wound contact layer
Facing:	Siliconized paper and plastic
Protective paper:	Cold seal paper coated with latex-based adhesive
Wrapping:	Carton

Shelf life: 5 years

### 3 LOGISTIC INFORMATION

	CU	SKU	Pallet
<b>GTIN:</b>	7310610064709	7310806470000	7312020064706
<b>Content:</b>	6 refills	10 CU	11 SKU/layer x 5 layers
<b>Size (w x l x h mm):</b>	125 x 100 x 90	130 x 513 x 190	Pallet 800 x 1200 x 1100
<b>Gross weight:</b>	205 g	2263 g	149 kg
<b>Packing material:</b>	Carton 25 g	Corrugated board 213 g	-

