

PRODUCT SPECIFICATION

REF 7221 CEDERROTH EYE WASH POCKET MODEL



1 GENERAL PRODUCT DESCRIPTION

Description

A personal eyewash bottle containing sterile, buffered, saline eye wash solution. The Cederroth Eye Wash Pocket Model is intended for washing away foreign particles and dangerous liquids, like alkali or acids, from an eye in an emergency first aid situation at work places.

The Cederroth buffered Eye Wash also has a neutralising effect on splashes from alkalis and acids. The effect is in general stronger for alkalis than for acids.

The Cederroth Eye Wash is intended for the first acute phase at the very accident site before getting to a plumbed or self-contained eye wash equipment or seeking medical attention.

The bottle is easily opened and ready to use by twisting off the dustcover. The integrated eye cup helps keeping the eye open during rinsing and ensures a generous flow without unnecessary waste. The bottle is shaped to fit into a pocket.



Accessory: Holster for Cederroth Eye Wash pocket model, REF 720300

Storage information

Storage in +2°C - +40°C

2 DETAILED PRODUCT INFORMATION

Content per consumer unit

<u>Product/Material</u>	<u>Qty</u>	<u>Dimensions (l x w x h mm)</u>
Cederroth Eye Wash Pocket Model Fill volume: 235ml	1	45 x 66 x 163

Material

Liquid	0.9% saline solution, borate buffered to neutral pH. Preservative free.
Bottle	Leak proof polypropylene bottle produced with Blow-Fill-Seal technique. Eye cup in transparent polypropylene. Dust cover in green hard plastics (ABS).

Shelf life: 4,5 years

3 LOGISTIC INFORMATION

	CU	SKU	Pallet
GTIN:	7310617221013	7310807221014	7312020012073
Content:	1 bottle	3 x 6= 18 CU	11 SKU/layer x 6 layers
Size (w x l x h mm):	45 x 66 x 160	157 x 468 x 173	Pallet 800 x 1200 x 1188
Gross weight:	294 g	5543 g	391 kg
Packing material:	-	Corrugated board 172 g Plastic 79 g	-

