

PRODUCT SPECIFICATION

REF 291400 CEDERROTH SMALL FIRST AID CABINET



1 GENERAL PRODUCT DESCRIPTION

Description

A metal cabinet with a small first aid panel attached to the front door and an inner storage compartment with two shelves, containing first aid product refill packages.

The Medical device products inside are CE-labelled.

Storage information

Dry storage in room temperature. For maximum protection, keep in original package until use.



2 DETAILED PRODUCT INFORMATION

Content per consumer unit

<u>Product/Material</u>	<u>Quantity</u>		<u>REF (refill)</u>
	<u>Instant use</u>	<u>Storage and refill</u>	
Metal cabinet	1	-	-
Small first aid panel (incl. plaster dispenser)	1	-	-
Salvequick plastic plasters (40 pcs)	1	6	6036
Salvequick textile plasters (45 pcs)	1	6	6444
Cederroth 4-in-1	1	1	1910
Cederroth 4-in-1 Mini	3	1	1911
Protection kit	-	1	2596
Tweezers	-	1	1875
Scissors	-	1	270002

Material

Cabinet:	Green metal
Outer door:	Polystyrene (PS)
First Aid products:	Refer to Product specifications for each product

Shelf life: 5 years from production date

3 LOGISTIC INFORMATION

	CU	SKU	Pallet
GTIN:	-	7310802914003	-
Content:	1 cabinet	1 CU	8 SKU/layer x 3 layers
Size (WxHxD mm):	260 x 440 x 170	280 x 480 x 190	800 x 1200 x 990
Gross weight:	4260 g	4800 g	140 kg
Packing material:	-	540 g	-



Version: 05
Issued: 2018-02-28