

PRODUCT SPECIFICATION

REF 3239 SALVEQUICK WOUND CLEANSER SACHET



1 GENERAL PRODUCT DESCRIPTION

Description

Sterile sachets made of non-woven tissues, which are saturated in a wound cleansing liquid. The tissues are packed in a well-sealed foil, which prevents leakage and guarantees sterility. Salvequick Wound Cleanser Sachet is intended for mechanical cleansing of minor wounds when clean tap water is not available.

The product is a medical device, class IIa.

Storage information

Dry storage in room temperature. For maximum protection, keep in original package until use.



2 DETAILED PRODUCT INFORMATION

Content per consumer unit

Product

Salvequick Wound Cleanser

Tissue

Sachet

Quantity

1000

Dimensions (mm)

100 x 140

60 x 58

Material

Tissue: Nonwoven made of viscose
Liquid: Physiological saline solution. Fill volume: 2,7 ml
Sachet: Paper/LDPE/Aluminium/Surlyn®

Sterilization: Gamma irradiation

Shelf life: 5 years from production date

3 LOGISTIC INFORMATION

	CU (transport package)		
GTIN:	7310803239006		
Content:	1000 sachets		
Size (w x l x h mm):	326 x 135 x 159		
Gross weight:	4450 g		
Packing material:	Carton 229 g, Plastic 42 g		

