

# PRODUCT SPECIFICATION

REF 6754CAP SALVEQUICK DETECTABLE ADHESIVE BANDAGES,  
FINGERTIP AND REGULAR, STERILE



## 1 GENERAL PRODUCT DESCRIPTION

### Description

Detectable, sterile adhesive bandages (plasters), of two sizes, medium and fingertip, packed in a cardboard refill intended for the Cederroth (Salvequick) dispenser for adhesive bandages. The sterile bandages are protected by a paper cap.

The Salvequick Detectable bandages are detectable by most common metal detectors and the blue color also makes them easy to detect visually, which makes these plasters suitable for restaurants and the food industry. The plasters can easily be pulled out with just one hand and applied on the wound.

The adhesive bandages consist of a backing, evenly spread with an adhesive coating and with an absorbent pad fixed centrally to the dressing. When the adhesive bandage is pulled out from the package, the lower cover will remain behind. Thus one half of the bandage can be applied with one hand and after removing the cap and the rest of the cover the plaster is ready for use.

### Storage information

Dry storage in room temperature. For maximum protection, keep in original package until use.



## 2 DETAILED PRODUCT INFORMATION

### Content per consumer unit

<u>Product/Material</u>	<u>Quantity (pcs)</u>	<u>Dimensions (mm)</u>	
Refill packages	6	30 x 88 x 100	
<u>Per refill package:</u>		<u>Dressing</u>	<u>Wound Pad</u>
Medium plasters	15	72 x 19	24 x 12
Fingertip	15	66 x 39	28 x 12

### Material

Film/backing:	Blue, soft and flexible polyurethane/polystyrene foil
Adhesive:	Acrylate
Wound pad:	Laminate of viscose, polypropylene, polyethylene with polyethylene wound contact layer and a metal detectable aluminium foil.
Facing:	Siliconized paper and plastic
Protective paper:	Cold seal paper coated with latex-based adhesive and a paper cap
Wrapping:	Carton
<b>Shelf life:</b>	5 years

## 3 LOGISTIC INFORMATION

	CU	SKU	Pallet
<b>GTIN:</b>	7310610067540	7310806754001	7312020067547
<b>Content:</b>	6 refills	10 CU	11 SKU/layer x 5 layers
<b>Size (w x l x h mm):</b>	125 x 100 x 90	130 x 513 x 190	Pallet 800 x 1200 x 1100
<b>Gross weight:</b>	185 g	2063 g	138 kg
<b>Packing material:</b>	Carton 25 g	Corrugated board 213 g	-



Version: 11  
Issued: 2017-06-26