



PRODUCT SPECIFICATION



Product: HOLSTER FOR
CEDERROTH EYE WASH POCKET MODEL

Reference code: 720300



1. General description

A green nylon pouch with a half-oval shaped covering flap. The pouch has a transparent window on the front side and a belt loop on the backside. The belt loop has both horizontal and vertical openings.

2. Intended use

The product is intended to be used for carrying a bottle of Cederroth Eye Wash Pocket Model e.g. in a belt.

3. Technical description

Product dimensions

(l x w x h) cm: 16 x 8 x 4,5

Material

Outer Fabric: Green Nylon 420D
Lining: Polyurethane
Transparent window: Transparent EVA

4. Composition

Consumer package

<u>Ref#</u>	<u>Product/Material</u>	<u>Pcs</u>
51618203	Holster in PE bag	1

Transport package

Number of consumer packages in transport package:	10
Dimensions (l x w x h) cm:	21 x 16 x 11
Weight (kg):	0,45

5. Labeling

On product/ consumer package

Cederroth address, telephone- and faxnumber, web-address
Product name, Reference code
Lot-number, EAN-code

On transport package

Cederroth address, telephone- and faxnumber, web-address
Product name, Reference code
Lot-number, EAN-code, Quantity

6. Storage conditions

The product should be stored in the original package for maximum protection until time of use.

Authorization: