

PRODUCT SPECIFICATION

REF 490750

SALVEQUICK PLASTER DISPENSER,
BLUE DETECTABLE



1 GENERAL PRODUCT DESCRIPTION

Description

A dispenser for adhesive bandages (plasters) which protects the plasters from dust, dirt and moisture until time of use. When a plaster is needed it can easily be pulled out from the dispenser and applied on the wound using one hand only. The plaster refills are locked in the dispenser and to remove them it is necessary to use the red key that is delivered with the product, a feature that works as an anti-theft system.



Plasters inside are medical device products and CE-labelled.

Storage information

Dry storage in room temperature. For maximum protection, keep in original package until use.

2 DETAILED PRODUCT INFORMATION

Content per consumer unit

<u>Product/Material</u>	<u>Quantity</u>	<u>REF (refill)</u>
Dispenser and lid	1	-
Blue detectable plaster refill	2 refill with 35 plasters	6735CAP
Key to dispenser	1	-

Material

Green parts:	Polystyrene and ABS-plastic
Transparent lid:	Polycarbonate
Sealing profiles:	EPDM rubber (synthetic rubber)
Blue detectable plasters:	Refer to Product specifications for REF 6735CAP

Shelf life: 5 years from production date

3 LOGISTIC INFORMATION

	CU	SKU	Pallet
GTIN:	7310614907507	7310804907508	7312024907504
Content:	1 dispenser	20 dispensers	4 SKU/layer x 4 layers
Size (WxHxD mm):	226 x 54 x 114	512 x 316 x 236	800 x 1200 x 1094
Gross weight:	258 g	5952 g	120 kg
Packing material:	N/A	Corrugated board 392 g Plastic tray 400g	-



Version: 02
Issued: 2018-02-28