PRODUCT SPECIFICATION

REF 490920 CEDERROTH FIRST AID STATION



1 GENERAL PRODUCT DESCRIPTION

Description

A First Aid station containing a variety of first aid products, intended to be used in first aid situations. All products are easily available and protected against dust and moisture by transparent lids. Space for additional products is provided behind a lid with first aid instructions. Adhesive bandages (plasters) and wound cleansers are placed in dispensers for easy access and theft-proof storage. A special key is provided for refill replacement.

Medical devices in the station are CE-labelled.

Storage information

Dry storage in room temperature. For maximum protection, keep in original package until use.



2 DETAILED PRODUCT INFORMATION

Content per consumer unit

| <u>Product</u> | <u>Qty</u> | <u>REF</u> | <u>Product</u> | <u>Qty</u> | <u>REF</u> |
|-------------------------------------|------------|------------|---------------------------------|------------|------------|
| | | (refill) | | | (refill) |
| First Aid station (empty) | 1 | - | Salvequick Wound Cleanser 20-p | 1 | 323700 |
| Cederroth 4-in-1 | 2 | 1910 | Protection kit | 1 | 2596 |
| Cederroth 4-in-1 Mini | 3 | 1911 | Cederroth Burn Gel Dressing 2-p | 1 | 901900 |
| Salvequick Plastic Plasters (45) in | 1 | 6036 | Cederroth Eye & Wound | 1 | 726000 |
| dispenser | | | Cleansing Spray | | |
| Salvequick Textile Plasters (40) in | 1 | 6444 | First Aid Instruction (label) | 1 | - |
| dispenser | | | Key to the dispensers | 1 | 51616620 |

Dimensions (w x d x h): (station) 558 x 286 x 117 mm

Weight: (full station) 2840 g

Material

First aid station: Green parts: Polystyrene, ABS-plastic

Transparent lid: Polycarbonate

Sealing profiles: EPDM rubber (synthetic rubber)

First aid products: Refer to Product specifications for each product

Shelf life: 5 years from production date for each product in the station.

3 LOGISTIC INFORMATION

| | CU | SKU | Pallet |
|----------------------|-----------------|-----------------|-------------------------|
| GTIN: | 7310804909205 | 7310804909205 | 7312024909201 |
| Content: | 1 station | 1 CU | 11 SKU/layer x 3 layers |
| Size (w x l x h mm): | 294 x 572 x 123 | 294 x 572 x 123 | 800 x 1200 x 1032 |
| Gross weight: | 3153 g | 3153 g | 129 kg |
| Packing material: | Carton 313 g | - | - |



Version: 03 Issued: 2018-03-01